

OCULAR SURGERY NEWS



The voice of 65,000 ophthalmologists worldwide™

Vol. 20, No. 12 • June 15, 2002

Two-year follow-up: radiosurgery better than laser

After 2 years of study, surgeons find that radiosurgery causes less heat damage than laser in blepharoplasty fat resection.

by David Bruce Welch, MD,
and Paul Bryar, MD

Special to Ocular Surgery News

An optimal instrument for the dissection, isolation, vascular coagulation and resection of fat pockets in blepharoplasty surgery must minimize the extent of excessive lateral heat diffusion. Too little heat diffusion during lipoplasty may result in vision-threatening hemorrhage. Too much lateral heat diffusion distorts tissues, interfering with precise and appropriate fat resection. Excessive lateral heat spread increases the degree and duration of postoperative inflammation, prolonging the recovery and increasing the risk of cicatrix formation.

Surgeons have used electrocautery, high-fluence carbon dioxide (CO₂) lasers, electrosurgery and ultrahigh radio frequency-modulated surgery (RFMS) for fat resection in blepharoplasty. Electrocautery uses electrical energy to heat a filament that becomes a red-hot electrode; incisions of human tissue are associated with extensive lat-

eral heat damage. Pulsed, high-fluence CO₂ lasers require pulse stacking with increasing of lateral heat diffusion and

“When tissue is placed between the electrodes, contact by the active electrode tip results in an incision.”

— David Bruce Welch, MD,
and Paul Bryar, MD

possible char formation via a cumulative heat sink phenomenon.

Electrosurgery and RFMS use fine delivery tips (active electrode) to transmit a current to a passive electrode (antenna). When tissue is placed between the electrodes, contact by the active electrode tip results in an incision. The incision occurs because of the tissue's natural resistance. A constant (non-modulated) energy waveform facilitates surgical cutting with enough thermal effect to close fine blood (less than 0.5 mm in diameter).

RFMS was first introduced in 1970. By using ultrahigh-frequency, low-temperature radio waves (4 MHz) a microfiber wire electrode can develop pressure-free incisions on tissue contact with no char to the incision margins and base. Larger-caliber blood vessels (greater than 0.5 mm in diameter) in both orbital and preaponeurotic

fat pockets must be coagulated prior to fat resection. This is necessary in both RFMS or CO₂ laser-assisted fat resection during cosmetic blepharoplasty.


Two-year follow-up

During the last two years we have used either RFMS (Surgitron DF Surge, Ellman International) or a pulsed CO₂ laser (UltraPulse, Lumenis) for fat resections on patients undergoing routine cosmetic transconjunctival lower eyelid


blepharoplasty.

Some yellowing of the incision margins and base were consistently observed with the laser incisions. Longer dwell time is necessary with the CO₂ laser to excise fat, compared with RFMS. We compared histopathologically the extent of lateral thermal damage in fat resected during lower eyelid blepharoplasty using these two surgical modalities.

Following the prolapse of redun-
(Radiosurgery results, continued on page 80)



1 Surgical margin of tissue excised with the CO₂ laser.



2 Surgical margin of the tissue excised with RFMS.

High-Frequency Radiosurgery

A New Patented Technology in Oculoplastic Surgery


ellman
Surgitron
4.0 Dual Radiofrequency
The Future of Surgery

- Endoscopic Brow Lifts
- Bipolar Hemostasis
- Punctal Occlusion
- Precise Incisions
- Frozen Sections
- Lesion Removal
- Bleph - Ptosis
- Trichiasis
- Pterygia
- DCR



OCULOPLASTIC TECHNIQUES USING HIGH FREQUENCY RADIOSURGERY

Video CD Demo By Jay Older, M.D.



The value of radiosurgery in oculoplastics: “The use of radiofrequency has markedly improved the performance of oculoplastic surgery”. With 4.0 MHz high frequency radiosurgery, blepharoplasties and other eyelid procedures can be done with minimal bleeding, minimal swelling, and maximum patient satisfaction. Dr. Older has found that this results in less surgical time and less post-op complications.



For more information:

(800)835-5355 • (516)569-1482 • www.ellman.com

ellman international, inc. 3333 Royal Avenue, Oceanside, NY 11572 • Fax: (516)569-0054

Radiosurgery results

(continued from page 79)

dant central preaponeurotic fat during lower eyelid blepharoplasty, the fat from one lower eyelid was excised using the CO₂ laser in a quasi-contin-

“Damage was determined by observing specific tissue changes such as cell lysis, irregularities in intact cell membranes and thickening and disruption of intercellular collagen.”

— **David Bruce Welch, MD, and Paul Bryar, MD**

uous wave mode at 8 W using a 0.2-mm focal tip. Fat from the contralateral lid was excised using RFMS in a cutting mode at 65 W with a fine-point needle electrode (Vari-tip).

After resection the fat was placed in 10% neutral buffered formalin and paraffin blocks were made. Then 5- μ m-thick sections were cut, placed on glass slides and stained with hematoxylin and eosin. The pathologist was masked regarding the instrumentation used for the resection of each specimen.

Results

Previous studies have described structural changes induced by various lasers on animal and human tissues. Nuclear and cytoplasmic changes have been observed. Nuclear changes include pyknosis, changes in chromatin patterns and nuclear necrosis. Cytoplasmic changes include loss of cell structures with preservation of the cell membrane. Coagulation necrosis and reduction of cells to an amorphous mass can also be seen.

Thermal damage at the margins of the resected tissue were assessed by measuring the depth of thermal injury. Damage was considered positive if any of the histologic changes described were present. This was determined by observing specific tissue changes such as cell lysis, irregularities in intact cell membranes and thickening and disruption of intercellular collagen. The amount of lateral tissue damage in each specimen was measured by a micrometer under 20X magnification power.

Two photomicrographs of the pieces of resected fat are shown. The

surgical margin is at the top of each of the photomicrographs. The tissue that was excised with the CO₂ laser (figure 1) shows more extensive thickening of

“Minimizing lateral thermal damage at the edges of fat resected during blepharoplasty should facilitate healing.”

— **David Bruce Welch, MD, and Paul Bryar, MD**

the collagenous septae and cell membranes than the tissue excised with RFMS (figure 2). These thermal changes are more prominent at the edge of the resection and extend deeper into the tissue when comparing tissue to that resected by RFMS. An average of 257 μ m of damage is observed with the CO₂ resected fat vs. an average of 148 μ m of damage with RFMS (H&E stain 10X magnification).

Ultrahigh RFMS resection of fat during lower eyelid blepharoplasty produced less thermal damage than the CO₂ laser. Prior to fat excision, medium-to-large blood vessels (greater than 0.5 mm in diameter) in harm's way must be selectively coagulated

prior to fat resection.

When using the CO₂ laser, these vessels were closed by using a defocused beam proximal to the incisional plane. When using ultrahigh RFMS, the vessels were clamped with steel forceps and were coagulated by conducting energy along the forceps by contact to the forceps with the Vari-tip.

Minimizing lateral thermal damage at the edges of fat resected during blepharoplasty should facilitate healing. The natural resistance of the fat to be resected by ultrahigh RFMS facilitates resection with less peripheral thermal damage than with the pulsed CO₂ laser.

OSN

For Your Information:

David Bruce Welch, MD, an assistant professor of clinical ophthalmology at Northwestern University Medical School and an attending oculoplastic surgeon at Northwestern Memorial Hospital, can be reached at 900 North Michigan Ave. (annex), Suite 1430, Chicago, IL 60611; (312) 255-0010; fax: (312) 337-5852. Dr. Welch receives honoraria from Ellman International. **Paul Bryar, MD**, an instructor of clinical ophthalmology at Northwestern University Medical School and an attending ophthalmologist, ocular pathologist at Northwestern Memorial Hospital, can be reached at 675 North St. Clair St., Suite 15-150, Chicago, IL 60611; (312) 695-8150; fax: (312) 695-3652. Dr. Bryar has no direct financial interest in the products mentioned in this article, nor is he a paid consultant for any companies mentioned.

■ Ellman International Inc. can be reached at 3333 Royal Avenue, Oceanside, NY 11572-3625, (516) 569-1482; (516) 594-3333; (800) 835-5355; Fax: (516) 569-0054 www.ellman.com

■ Lumenis can be reached at 2400 Condensa St., Santa Clara, CA 95051; (408) 764-3000; fax: (408) 764-3660.

Know Me?

Patients Do

CRSQA is a patient/consumer organization committed to help refractive surgery patients find the best available doctors. That is one of the reasons I decided to become a CRSQA Certified Refractive Surgeon. Patients rely on CRSQA to help them choose a surgeon. I rely on CRSQA's certification to help inform patients of my qualifications. Peyton Neatrour, MD - Virginia Beach, VA

Council for Refractive Surgery Quality Assurance
Nonprofit Patient Education and Refractive Surgeon Certification
www.usaeyes.org
800/USA-EYES



Register for the OSN SuperSite



On the Web with OSNSuperSite.com

The OSN SUPERSITE consolidates volumes of ophthalmic information into one manageable resource. Finally, an online resource that puts the power of the Web to work for you.

Here's how to register:

1. Call up the OSN SUPERSITE home page at: **OSNSuperSite.com**
2. When attempting to view the full text of an article, first-time users will be automatically prompted to Register on the OSN SUPERSITE Web site.
3. Click on Register under Welcome New Users.
4. Fill out the standard Registration Information. *Note:* the e-mail field must be filled out.
5. If you're an ophthalmologist, tell us your primary practice focus.
6. Your Username and Password will be stored in a cookie on your computer. Please be sure cookies are enabled on your browser. If you would like to access the site from a different computer, you do not need to re-register. Simply login with your name and password. If you cannot remember your password, you may retrieve it by providing your name and e-mail address.



OSNSuperSite.com